

## Product Description


Device used for providing a specialized source of light for illumination of a site of medical intervention. It provides a high intensity, high colour rendering field of light that minimizes shadows and the emission of heat. Adjustable Operating field diameter with focus



Automatic switch from mains to batteries in case of power failure  
With battery status indicator  
Removable autoclavable handle  
Power supply: The supplier has to supply the right voltage, plug and socket type that fit the country it is delivering the goods to. The country of destination would be stated in the UNFPA's



Mobile on heavy sturdy stand,  
On 4 antistatic ball-bearing swivel castors, 2 with breaks  
Spring balanced articulating arm  
Low centre of gravity for optimal stability and reach



## Main Technical Data

YD01-5LED	Shadowless Operating lamp
Illuminance	≥100,000
Colour temperature	4800±300K
Import power	150VA
Illuminance depth	≥500mm
Total irradiance	547W/ m²
Size of light field	160-220mm
Lamp head diameter	500mm
Rated power of bulb	3.2V/1W(32pc led bulb )
Power supply voltage	220V/50Hz
Bulb	cold light